

CLINIS REAGENTS STORAGE

ASSAY	REAGENT	TEMPERATURE [°C]	AFTER OPENING
Total cholesterol	cuvette	15-30	4M
	enzyme	2-8	-
Glucose	cuvette	15-30	4M
	enzyme	2-8	-
Triglycerides	cuvette	15-30	6M
	enzyme	2-8	-
HDL	cuvette	2-8	6M
	enzyme	2-8	-
Hemoglobin	cuvette	15-30	6M
Haematocrit	cuvette	15-30	6M
Erythrocytes	cuvette	15-30	6M
BITEST	cuvette	15-30	4M
	Enzyme 1	2-8	-
	Enzyme 2	2-8	-
Uric acid	cuvette	15-30	6M
	enzyme	2-8	-
FORT	R1	15-30	-
	R2	15-30	-
FORD	C1	15-30	-
	S1	15-30	-
	S2	15-30	-
	S3	15-30	-
HbA1c	R1	2-8	-
	R2	2-8	-
	R3	2-8	-
ALT/AST	cuvette	2-8	-
	Reagent R1	2-8	-
	Reagent R2	2-8	-

Callegari S.r.l. is a Catellani Group Company

Headquarter – Castello dei Diritti, Via Eugenio Ugolotti, 1 43122 Parma/ Manufacturing Site – Via Luigi Natale Vernazzi, 13/A 43122 Parma

VAT No IT01615440342 - REA No PR-166829 – Share Capital € 516.000,00 Fully Paid